



Environmental Control Systems

User manual for IR-transmitters

Pico4, Pico8 and Pico16



Revisions

Changes in document will be listed here.

Revision list

Date	Name	Description
15.11.2015	Bent Sollid	First version. Based on Norwegian version 01.12.15.
06.06.2016	Bent Sollid	Ch 3.3.1: added details for Android.
11.11.16	Bent Sollid	Ch 3.3.2: added details for Apple ref iOS 10.0.1.
01.05.17	Bent Sollid	Ch 3.: figure exchanged.
01.10.19	Bent Sollid	Ch 2.7 and 8.3: details about battery and charging.

Contents – IR-transmitters Pico4, Pico8 and Pico16

1	Introduction.....	4
1.1	Pico4, Pico8 and Pico16 in general	4
2	User manual.....	4
2.1	Light indicators and connections	4
2.2	Ordinary use	5
2.3	Macro.....	5
2.4	Page 2 enabled	6
2.4.1	Maximum no of channels on Pico4 and Pico8 when page 2 is enabled.....	6
2.5	Lock code key	6
2.5.1	Common code	6
2.6	Scanning	6
2.7	Battery	7
2.7.1	Empty battery	7
2.8	Outdoor use.....	7
3	Programming – on screen	8
3.1	Current consumption in programming-mode	8
3.2	Connections	8
3.3	Start-up	9
3.3.1	Android.....	9
3.3.2	Apple	9
3.3.3	Windows PC.....	10
3.4	What's on the screen.....	10
3.5	Explaining some settings	11
3.6	Copying IR-signals into the transmitter.....	13
3.6.1	Details about IR-light	14
3.7	Macro.....	14
3.8	Factory setting	14
3.8.1	Pico4.....	14
3.8.2	Pico8.....	15

3.8.3	Pico16.....	15
3.9	Backup.....	16
4	Programming – on IR-transmitters.....	16
4.1	Programming mode – enter and exit.....	16
4.2	Enable/disable common code	16
4.3	Copying IR-signals	16
5	Pico4 and Pico8 – keypad exchange	17
5.1	How to do it	18
6	Maintenance	18
6.1	Cleaning.....	18
7	Troubleshooting	19
7.1	If errors occur	19
8	Technical information	20
8.1	Battery	20
8.2	External switch	20
8.3	Recycling	20
8.4	Technical data.....	20
8.5	CE-markings	21
9	Appendix – programming map.....	22

1 Introduction

This user manual is meant for Picomed's IR-transmitters Pico4, Pico8 and Pico16. These handheld IR-transmitters have four, eight or eighteen keys and two pages. Pico4, Pico8 and Pico16 IR-transmitters may be programmed via a browser when i.e. a computer is connected to the transmitter with a USB cable. A few functions may be programmed by the use of the keys on the transmitters.

1.1 *Pico4, Pico8 and Pico16 in general*

The IR-transmitters are meant used as handheld transmitters and can be hung in a lanyard or mounted on a table or wheelchair. The integrated finger guide makes it easier for people with spasms to operate it.

- Pico4 and Pico8 has the same size but with different number of buttons.
- Pico16 has 18 keys: 16 for IR-transmitting, 1 key for page select and 1 for lock code.

Features:

- Learnable IR.
- Built in IR-signals for environmental control systems.
- 1 switch input for scanning.
- Macro.
- Integrated finger guide.
- Continuous or single IR-transmission when press & hold on key.
- Programmable: sound on keypress, spasm filter, parameters for scanning, backup.

Pico4, Pico8 and Pico16 IR-transmitters can learn many different IR-signals on the market, and can easily be programmed to be used on different remote controlled systems. It is operated with soft and low force keys. Light indicators are used for indicating IR-transmission and when programming the unit.

2 User manual

Content on delivery:

- IR-transmitter.
- Battery charger and USB-cable.
- USB cable «OTG» for interconnection between IR-transmitter and unit for programming.
- Lanyard.
- User manual.
- Pico4 and Pico8: extra keypad with arrows: to replace original keypad if customer want it instead of original keys with numbers.

Optionally, an external switch with 3,5 mm mono jack can be ordered for scanning-mode.

2.1 *Light indicators and connections*

The figure below shows light indicators and connections on a Pico4. Pico8 and Pico16 has the same functions.



2.2 Ordinary use

Pico4, Pico8 and Pico16 are made as handheld transmitters and may be hung in a cord around the neck. Infrared (IR) light is transmitted from the front of the transmitter. If possible, the IR-transmitter should point towards the receiver when it is used. This will increase the range and make it easier for the user to see the light indicators. Obstacles between the transmitter and the receiver may stop IR-signals.



To operate the transmitter, press a key.

- A green light indicates IR-transmission.
 - If continuous IR-transmission is enabled, IR-signal will be repeated as long as key is pressed.
 - If continuous IR-transmission is disabled, IR-signal will be sent only once when key is pressed.
- If no IR is programmed on the key, a red light indicator will indicate it.
- If battery is low, a red flashing light will indicate it.
- If beeper is enabled, it will give a beep for each keystroke.
- If page 2 is enabled, the yellow page-indicator will illuminate some seconds at each keystroke. Yellow light on the left side is for page 1 while yellow light on the right side indicates page 2. More information about this in chapter 2.4.

A lanyard and a rubber-cover is included. To avoid dust and fluid spill inside the controller, it is recommended to attach the cover when connections are not in use.

2.3 Macro

A function named macro will be useful if there is a need of transmitting more signals with only one keystroke on the IR-remote or on an external switch connected to the remote control.

- A green light indicates IR-transmission.
- If battery is low, a red flashing light will indicate it when a key is pressed.
- If beeper is enabled, it will give a beep for each keystroke.

- A macro can be interrupted by pressing a key on remote control or external key connected.

2.4 Page 2 enabled

When page 2 is enabled, the page-select key will switch between page 1 and page 2. When not in use, the transmitter will remember the last page used. A yellow light will illuminate some seconds after each keystroke to indicate which page is selected: left is page 1 and right side is page 2.

- On the Pico4 and Pico8, the lower left key is the page-select key when page 2 is enabled. This key will not transmit IR when page 2 is enabled.
- On the Pico16 the lower left key, marked with two arrows, is dedicated as page select key. If page 2 is *disabled*, this key will not have any function.

2.4.1 Maximum no of channels on Pico4 and Pico8 when page 2 is enabled

- Pico4 has 3 keys on each page for IR-transmission, will be able to transmit 6 channels.
- Pico8 has 7 keys on each page for IR-transmission, will be able to transmit 14 channels.
- If lock-code key is programmed as page-independent (same IR-signal on page 1 and page 2), it will transmit the same signal on page 1 and page 2.
Pico4 will have a maximum of 5 channels.
Pico8 will have a maximum of 13 channels.

2.5 Lock code key

The lock code key transmits IR-signals in the same way as the other keys. Additionally:

- It may be programmed to be page independent. In this case, the same IR-signal will be transmitted regardless page 1 or page 2 is selected.
- It is able to transmit a «lock code», i.e. a more advanced code which is more secure than ordinary code. Normally, these codes is used for entrance doors.
- It is able to send a «common code». For more information, please refer to chapter 2.5.1.

2.5.1 Common code

The lock code key can be programmed to transmit a «common code» together with the transmission of a lock code. This signal is mainly used on common doors which is accessed by many people. With common code enabled, it is possible to use the same key for opening a common door (activated by common code from transmitter) and for entrance doors (activated by lock code from transmitter). The IR-transmitter transmits both common- and lock code when pressing the key once. No special action necessary to choose one or another signal.

2.6 Scanning

The scanning mode may be used when a person is not able to use the ordinary keys. The IR-transmitters may be controlled by one external switch connected with a 3,5 mm mono jack. It is also possible to use ordinary keys when an external key is connected.

- If an ordinary key is used when scanning is activated, the scanning sequence is cancelled. Picomed's transmitter will respond to this as an ordinary keystroke.
- Press & hold on external switch will not be cancelled by pressing an ordinary key.

- Press & hold on ordinary key will not be cancelled by pressing external key.
- Keys without a programmed IR-signal is not included in scanning sequence.

When pressing an external switch connected to the IR-transmitter, it will start the scanning sequence.

- Scanning starts on first key on last page used.
- A light in each key will illuminate for a chosen time.
- When all keys on a chosen page is scanned through, it will start over again with the first key on the same page.
- When external switch is pressed again, the IR-transmitter will transmit the IR-signal belonging to the key that is illuminated at the moment.
- The green and red lights will light up as described in chapter 2.2.
- When IR-transmission is finished, scanning sequence will halt on actual key for 3 seconds before it proceeds. This makes it easy to activate new IR-transmissions on same key.
- If page 2 is enabled, page select key will be included in the scanning sequence.
- Scanning stops after a pre-set number of scanning's through all active keys. The pre-set number of scanning's may be changed in programming mode.

If beeper is enabled in scanning mode:

- A long beep for each keystroke on external key.
- A short beep when going from one key to the next key.
- A long beep when going from last key to the first key on a page.

2.7 Battery

Pico4, Pico8 and Pico16 has a rechargeable battery.

- Red flashing for 2 seconds when transmitting IR indicates low battery.
- Recharge battery with the included battery charger.
- When charging: green slow flashing.
- When fully charged: steady green light.
- Typical time of charging from empty battery is 4-5 hours.
- Charging to be done in room temperature. Minimum temperature shall be 10 °C.
- The transmitter can be used when charger is connected. In this case, the green charging indication will be disabled some seconds when it transmits IR. Information about light indicators may be found in chapter 2.2.
- A cold battery will normally have lower performance than a room tempered battery.

2.7.1 Empty battery

The battery shall not be empty when stored. Recharge battery before the IR-transmitter is placed for storage for a long time. If the IR-remote is not used for a long time, it is recommended to recharge it once a year.

No settings, including IR-signals, will be deleted if the battery is empty.

2.8 Outdoor use

Picomed IR remote controllers can be used outdoors, but they should be protected against water and fluids. Avoid spilling of drinks since this can result in damage to the remote

controller. In a case of fluid entering the remote controller, the controller should be removed from the fluid at once, and any liquid drained out of the case. Place the remote controller at room temperature in a dry and airy place with the cover removed. Wait at least 24 hours and check if it still works.

The IR remote controller will work and transmit IR signals also in severe cold, but the use is limited by the fact that the battery will have low efficiency at low temperatures. This may result in a low battery indication that disappear when it is back in room temperature.

3 Programming – on screen

The IR-transmitters Pico4, Pico8 and Pico16 can be programmed by the use of an external unit with a screen. I.e. a smartphone, tablet or PC with a web-browser. With a USB cable connected to the transmitter, a webpage stored inside the transmitter can be shown in the browser. No data will be downloaded from internet, but some units needs an internet connection available to be able to load a web page via the USB cable.

Some functions in the IR-transmitters may be programmed directly on the transmitter without the use of an extra unit with screen. More information about this is found in chapter 4.

3.1 Current consumption in programming-mode

- It is recommended to recharge both the transmitter and smartphones/tablets before programming starts.
- It is recommended to disconnect USB cables from the transmitter when programming is finished.
- If low battery occurs while programming, it will be indicated by a red flashing light every minute. In this case, the battery has to be recharged.

3.2 Connections

There are two ways to connect the cable(s) depending on which unit to use when programming the IR-transmitters.

- The USB-cable supplied with the IR-transmitter is used for charging its battery and when connecting to a device with screen.
- This will not fit in some units, for instance Apple products. In such case, one have to use the charging cable supplied with your unit.



Charging cable




USB adapter cable
«OTG»

Normally, see figure below:

- Windows PC: use the charging cable supplied with your IR-remote control.
- Smartphone/tablet: connect the USB adapter cable (enclosed «OTG») to the transmitter Pico4-8-16.

- Connect your units charging/USB cable (or the one supplied with Pico4-8-16 if it fits) to the adapter cable.
- Some units need more external power from the USB cable than the IR-transmitter is able to supply. In such cases it may help to use an “Y-cable” connected to the charger. Contact your supplier for more info.



3.3 Start-up

3.3.1 Android

Android products use an app which has to be downloaded from Google-play. Search for «Picomed» and install the app «Picomed Starter».

- Connect cables as described in chapter 3.2.
- On some models, a USB icon will be shown on screen, then open the Picomed Starter app.
- Choose «Enable Connection» and enable «USB-internet sharing».
- Go one step backward in the app on your unit.
- Choose «Configure Product», the app will open the web-browser* on your device and download a web page from Pico4-8-16 and show it on your screen.
- See chapter 3.4.

*) If you get the message “Product found” but the web page is not opened:

- Some mobilephones want for more current than the remote control allows. In such case a “Y-cable” can supply the mobilephone from a charger at the same time as the remote is connected. Please contact manufacturer for more information.
- Some mobilephones have an antivirus program/app installed which refuse the web-page to be loaded via USB. You may have to deactivate or uninstall this app.
- Some mobilephones need to have “Roaming” activated before it is possible to load a web-page. You may have to activate mobile data before roaming can be activated.
- Some mobilephones will deactivate roaming when it connects to a wireless network (WiFi). In some cases it is possible to activate roaming after WiFi is connected.

3.3.2 Apple

Apple products will work without any installation on iPhone 4 or newer, and for iPad 3 or newer which has an integrated SIM-card with mobile data enabled. Tested on iOS 7, 8, 9 and 10.

Note: Apple has different behaviour on the different iOS versions. It will therefore be differences between the description below and the behaviour on iPhone on some iOS versions. In some cases the iPhone will ask twice for a “Ok” press. You will in some cases have to turn personal hotspot off and on twice depending on iOS version. One should wait a few seconds between each point in the list to let it have time for verifying connection etc.

- You have to turn on mobile data on the unit to get access to the “Personal hotspot” option. This option has to be enabled to let it use the USB cable for downloading an internet page from the Picomed transmitters.
- Enable mobile data: tap Settings – Mobile data – Mobile data = on.
- Turn off hotspot: tap Settings – Personal hotspot – Personal hotspot = off.
If this alternative is not available, it may be that your unit or your mobile supplier not has this option.
- Connect your device and the Picomed IR-transmitter with the USB cable.
- Wait for a message “Trust this computer” and tap “Trust”.
- Turn on hotspot: tap Settings – Personal hotspot – Personal hotspot = on.
If there is a question about which type, choose “Share only with USB”.
- Personal hotspot = off – wait a few seconds – Personal hotspot = on
- If there is a question about which type, choose “Share only with USB”.
- Wait for message on top of screen, “Personal hotspot: 1 connection” (blinking slowly).
- Return to main menu and open the internet browser in your product, enter address in browser: `picomed.local`
- See chapter 3.4.

3.3.3 Windows PC

A Windows PC has to install a driver the first time you are going to use it. When the computer is connected to internet, a driver will automatically be downloaded from Windows Update. This will normally take a minute.

- Connect USB cable between computer and IR-remote, open a web-browser and go to: `picomed.local`
- See chapter 3.4.

3.4 What's on the screen

Your unit will now load a web page from Picomed's transmitter. There is no need of internet data traffic to do this, since data is loaded via the USB cable from the transmitter. Please refer to chapter 9 to find an overview of the menus. The figure below shows the principle of operation.





Figure 1 Main page in the programming menu

By choosing «Keys on page 1», see figure above, you will get access to a menu to make settings for keys regarding page 1, and so on for the other settings available.

- To go back one level: swipe to the left or use the browser's «back» function.
- Some settings will ask for a confirmation if you change a value, see figure below, press yes/no to confirm or not.



- When changes are done, the transmitter confirm it by a low beep.
- Sliding switches are used to turn functions on/off.
 - Slide switch to the right to turn on a function.
 - Slide switch to the left to turn off a function.
 - Some sliding switches may control more than one value. If the sliding switch is positioned in the middle (neither on nor off), it means that some of these functions are on and some are off. If you slide it on or off, all functions controlled by this switch will be set either to on or to off.

3.5 Explaining some settings

Settings in brackets are factory settings. It is a graphical menu system, which in many ways are self-explanatory.


Topic

Keys on page 1/2	<p>This is access to settings for each page; choose which page to make settings on.</p> <ul style="list-style-type: none"> • Page 1 will always be enabled. • Page 2 can be disabled. On a Pico16, the page-select key will have no function. • Page 2 can be enabled. On Pico4 and Pico8 the lower left key will be a page select key without any IR function.
General settings	General settings for the IR-transmitter.
IR signal	Choose which IR-signals to transmit. It is possible to choose between built in IR-signals and to copy IR-signals from other transmitters. In chapter 3.8 the IR-signals which is preset in factory setting is listed.
Picomed	Picomed's IR-format («IR-language»).
User Code	Choose one of 32 user-codes, values from 00-31.
Channel	Choose one of 16 channels, values from 00-15.
Pre-set door phone	Choose a pre-set user code (05) and channels (13-15) used to control Picomed's door phone Picolo.
Pre-set phone	Choose a pre-set user code (05) and channels (00-15) used to control Picomed's telephone PicoPhone.
Gewa	Gewas IR-format («IR-language»).
Channel	Choose one of 128 channels, values from 00-127.
Learn/copy IR	Menu for learn/copy IR-signals from another IR-transmitter. Please also refer to chapters 3.6 and 4.3.
Test	Transmits the IR-signal, for testing in learning mode.
Continuous IR	If enabled, it will repeat the IR-signal until key is released. If disabled, it will transmit the IR-signal once even though a key is pressed for a long time.
Delete IR	Delete all IR for a key/page.
Individual key	Settings for an individual key on chosen page.
Choose key	Choose an individual key on chosen page.
Lock Code	If enabled, it will transmit lock-code. More information can be found in chapter 2.5. Only available for individual key «lock-code» placed in lower right corner on transmitter.
Common code before LC?	<p>If enabled, it will transmit a common-code before the lock code when lock code key is pressed. See chapter 2.5 and 2.5.1. Only available for individual key «lock-code» placed in lower right corner on transmitter.</p> <ul style="list-style-type: none"> • If key is programmed with a Picomed IR-signal, it will transmit Picomed common code. • If key is programmed with a Gewa IR-signal, it will transmit Gewa common code.

	<ul style="list-style-type: none"> If key is programmed with a copied IR-signal or if no IR is programmed for this key, it will transmit a common code depending of which country the product is configured for (see programming menu, general settings).
Picomed LC Gewa LC	Picomed's IR-format for lock code (LC = Lock Code). Gewas IR-format for lock code.
Lock key page2=page1	If enabled, the lock code key is page independent (off).
Macro	Make a macro. Functions from 1 to 10 keys can be executed in sequence when pressing one key.
Sound at keypress	Enable a short beep when a key is pressed (off).
Spasm filter Time setting	Enable spasm filter when a key is pressed (off). Spasm filter delay in seconds when key is pressed (1,5). Key <i>release</i> delay maximum value is 0,5 sec.
Scanning Sound when scanning Interval Number of rounds	Enable scanning when external key is connected (on). Enable sound when scanning is performed (off). Scanning time delay in seconds for each key (1,0). No of repetitions (2).
Factory setting	Reset all settings to factory values including deleting IR-signals copied into the transmitter, details in chapter 3.8.
Product status	Battery voltage. Note: when connected to a PC, the voltage will be higher due to charging supplied by PC.
Backup/restore	Enables backup to/from the transmitter, details in chapter 3.9.
Product Firmware	Information about versions. Uploading new firmware.

3.6 Copying IR-signals into the transmitter

- Each IR-signal has to be copied twice into Picomed transmitters.
- IR-signals programmed to a key will be deleted when a new signal is copied into it.
- Existing IR-signal will be deleted when the IR-remote has received a new IR-signal.
- IR-signals «in the room», i.e. sunlight, fluorescent lamps etc. will also be read into the IR-remote when copying IR-signals. The IR-receiver in the front of the Picomed transmitter should be shielded from other IR-signals than the ones from the IR-transmitter to be copied.
- If it is hard to make a useable copy of a signal, try to change distance (normally 0-10 cm) and angle between the IR-remotes.
- It is strongly recommended to press key short and concisely on the other remote control to avoid repeated signals due to a press & hold condition.



3.6.1 Details about IR-light

Infrared light is a not visible light that are in our environment all the day. Both as a natural radiation from the sun and as human made radiation from i.e. electrical light. When copying IR-signals to Pico4, Pico8 or Pico16, the IR-remote controllers are not able to separate the IR-light from the sun or electrical lights from the IR-signals from another IR-remote controller. Therefore, these IR-signals also will be copied into the Picomed controller. To avoid this situation it is recommended to shield the area over and between the Picomed IR-transmitter and the original IR-remote controller from other IR-signals while copying signals. This can be done by covering the area over the two of them with your hand, a paper, a book etc.

If not, the IR-receiver receives the IR-light in the room and stores it including unwanted signals. It may therefore finish the IR-training process based on unwanted signals before the IR-transmit key is pressed on the other remote controller.

3.7 Macro

The macro function makes it possible to transmit several IR-signals with only one keystroke.

- To be programmed from web-browser.
- Macros may be programmed on both pages, even though page 2 is disabled.
- Page select key on page 2 cannot be configured with a macro.
- There is no limit of how many keys which may be configured with macros, although if many keys are used for macro, fewer keys will be available to be used in macros.
- Each macro consists of a number of commands with info about which page and which key IR shall be sent from.
- A macro can have a maximum of 10 commands.
- All keys may be used in a macro.

3.8 Factory setting

When a factory setting is done, it will:

- Delete all IR-signals.
- Program IR-signals as shown in tables below.
- Program factorysettings for all parameters as shown in parantheses in chapter 3.5. «NOR» and «SVE» in table refer to Norwegian or Swedish configuration in the IR-remote. Configuration of your unit can be found in the general settings menu – Product Firmware.
- In the tables, LC = LockCode.

3.8.1 Pico4

- If page 2 is enabled, then key 3 will act as a page select key. If disabled, key 3 can transmit IR-signals in the same way as the other keys.
- A lock code can be programmed on Key 4.
- Key 4 may be programmed to transmit the same IR-signal on page 1 and 2; the key will be “page independent”.

Key	NOR page 1	NOR page 2	SVE page 1	SVE page 2
1	Picomed 5-13	< Empty >	Gewa 16	<Empty>
2	Picomed 5-14	< Empty >	Gewa 17	<Empty>
3	Picomed 1-1 ¹	Page select key	Gewa 18 ¹	Page select key
4	Random Picomed LC	<Empty>	Random GLC	<Empty>

Note 1: key 3 will act as a Page select key on page 1 and 2 if page 2 is enabled.

3.8.2 Pico8

- If page 2 is enabled, then key 7 will act as a page select key. If disabled, key 7 can transmit IR-signals in the same way as the other keys.
- A lock code can be programmed on Key 8.
- Key 8 may be programmed to transmit the same IR-signal on page 1 and 2; the key will be “page independent”.

Key	NOR page 1	NOR page 2	SVE page 1	SVE page 2
1	Picomed 5-13	<Empty>	Gewa 16	<Empty>
2	Picomed 5-14	<Empty>	Gewa 17	<Empty>
3	Picomed 1-1	<Empty>	Gewa 18	<Empty>
4	<Empty>	<Empty>	Gewa 19	<Empty>
5	<Empty>	<Empty>	Gewa 20	<Empty>
6	<Empty>	<Empty>	Gewa 21	<Empty>
7	<Empty> ¹	Page select key	Gewa 22 ¹	Page select key
8	Random Picomed LC	<Empty>	Random GLC	<Empty>

Note 1: key 7 will act as a Page select key on page 1 and 2 if page 2 is enabled.

3.8.3 Pico16

- The lower left key is a dedicated page select key.
It is not possible to transmit IR-signals with this key.
If page 2 is disabled, this key will have no function.
- The lower right key is a dedicated lock code key and will always be page independent.

Key	NOR page 1	NOR page 2	SVE page 1	SVE page 2
0	Picomed 5-0	<Empty>	Gewa 25	<Empty>
1	Picomed 5-1	<Empty>	Gewa 16	<Empty>
2	Picomed 5-2	<Empty>	Gewa 17	<Empty>
3	Picomed 5-3	<Empty>	Gewa 18	<Empty>
4	Picomed 5-4	<Empty>	Gewa 19	<Empty>
5	Picomed 5-5	<Empty>	Gewa 20	<Empty>
6	Picomed 5-6	<Empty>	Gewa 21	<Empty>
7	Picomed 5-7	<Empty>	Gewa 22	<Empty>
8	Picomed 5-8	<Empty>	Gewa 23	<Empty>
9	Picomed 5-9	<Empty>	Gewa 24	<Empty>
10/A/*	Picomed 5-10	<Empty>	Gewa 26	<Empty>
11/B/#	Picomed 5-11	<Empty>	Gewa 27	<Empty>
12/C	Picomed 5-12	<Empty>	Gewa 28	<Empty>
13/D	Picomed 5-13	<Empty>	Gewa 29	<Empty>
14/E	Picomed 5-14	<Empty>	Gewa 30	<Empty>
15/F	Picomed 5-15	<Empty>	Gewa 31	<Empty>

S	Page select key	Page select key	Page select key	Page select key
L	Random Picomed LC	¹	Random GLC	¹

Note 1: lower right key is page independent (lock code key).

3.9 Backup

- It is possible* to store all settings from the Pico4, Pico8 and Pico16 in a data file.
- A backup-file may be restored to the Pico4, Pico8 and Pico16.
- It is only possible to restore a backup to a similar IR-transmitter as it was made from.
- In the backup menu, there is a Compatibility ID that has to be equal in both transmitters if a backup/restore operation shall be possible.

*) Due to Apple's handling of files and folders this function is not available from Apple units.

4 Programming – on IR-transmitters

Some programming may be done directly on the IR-remote without the use of an extra unit with a screen. Beeper will be enabled when this programming mode is active.

4.1 Programming mode – enter and exit

- Enter programming mode: press & hold keys 1, 2 and 3 simultaneously for 3 seconds.
- Exit programming mode: press & hold key «1» for 3 seconds.
- The beeper gives a long beep to indicate that programming mode is entered.
- This programming mode will always start on page 1 indicated by left light indicator flashing.
- If low battery, the red light indicator will flash for 2 seconds and then the transmitter will leave programming mode.
- The transmitter will leave programming mode after 60 seconds without any keystrokes.
- Light indicators will be turned off and beeper will give a short beep when leaving programming mode.

4.2 Enable/disable common code

In programming mode, it is possible to enable/disable transmission of common code. Press & hold the lock-code key for 3 seconds to change setting. For information about common code, please refer to chapter 2.5.1.

- If disabled, this is indicated by a short beep, a blink in the light in the lock-code key and a blink in the green light indicator.
- If enabled, this is indicated by a long beep, a long blink in the light in the lock-code key and a long blink in the green light indicator.
- Which common code to be transmitted is dependent of the configuration of the product.
If set to "Norway": IR = Picomed, user code 1, channel 1 (earlier designated NA1-1).
If set to "Sweden": IR = Gewa, channel 16.

4.3 Copying IR-signals


In this programming mode, IR may be captured from other IR-remote controls. Also refer to chapter 3.6 and 3.6.1.

- The yellow light-indicator indicates which page is chosen.
- To switch page, press & hold page select key (lower left) for 3 seconds.

- IR capture can be done on page 2 even if page 2 is disabled.
- IR capture can be done on page select key on page 1 even though page 2 is enabled.
If page 2 is enabled, the lower left key will not transmit any IR because it is used as a page select key.
- IR capture cannot be done on page select key on page 2.

How to...

- Each IR signal has to be captured twice.
- Position Pico4-8-16 and the other IR-transmitter facing each other. It should be about 2-3 cm between the two of them.



- Press shortly on the key on Pico4-8-16 to be programmed. The key will lit, green light indicator will lit, yellow page indicator light will flash slowly.
- Press shortly on the key to capture IR from on the other remote controller. When first IR signal is captured there will be a short delay, then the green and yellow indicator lights will flash synchronously.
- Press shortly on the key on the other remote controller to capture IR from.
When first IR signal is captured, green and yellow indicator lights will flash synchronously.
- Press shortly once again on the key on the other remote controller to capture IR from.
When the second IR signal is captured, light in the key will be turned off, green light indicator will flash 3 times, and then the yellow light indicator continue flashing slowly.

Repeat sequence for other keys as needed.

If the Picomed controller does not receive IR signals within 20 seconds, it will exit to programming mode. This is indicated by the yellow flashing light indicator. Leave programming mode as described in chapter 4.1.

5 Pico4 and Pico8 – keypad exchange


Pico4 and Pico8 have alternative keypads.

- Pico4 (included): keypad with 4 arrows, may be rotated to point up-down or left-right.
- Pico8 (included): keypad with arrows, numbers and key-page select.
- Pico8 (not included): can combine to keypads from Pico4, arrows and/or numbers.



5.1 How to do it

- Turn remote controller up-side down, use a Torx tool size T10 and remove two screws.
- Keep front- and back part together, turn transmitter back with front side up.
- Remove front cover including keypad. Avoid touching electronic card.
- Replace keypad with the new one.
Pico4: keypad may be rotated if needed, see figure.
Please note that keypad shall fit into guiding marked on figure with a ring.
- Assemble in opposite order as disassembled and test functionality.



6 Maintenance

Picomed's IR transmitters Pico4, Pico8 and Pico16 are battery operated and shall be charged with the battery charger included. Normally there is no need of regularly maintenance. There are no fuses etc. to replace. It has to be cleaned when necessary.

6.1 Cleaning

Cleaning of the remote controller should only be done with a moist cloth. It may be used a mild detergent. It must not become wet or exposed to excess fluids.

7 Troubleshooting

In the case of a failure, the following table should be checked.

Symptom	Possible error	Correction
Transmitter does not respond	Empty battery.	Connect battery charger.
Transmitter has a short range or does only work when pointing directly against the receiver.	Low battery. Wrong placement of receiver. An obstacle is blocking the IR-signal.	Connect battery charger. Check if the receiver has a position where it can receive IR-signals without being disturbed. See above. Look for obstacle between transmitter and receiver.
Transmitter does work for a short time on some or all channels.	Transmitter not correct programmed – programmed without continuous IR.	Reprogram: enable continuous IR.
Transmitter does only function on some keys.	Transmitter not correct programmed – programmed with wrong codes.	Reprogram. Check if green light indicator illuminates when key is pressed.
When IR-capture: transmitter exit IR-capture before IR-signal is sent.	IR-noise from environment.	Avoid IR-signals from bulbs, sun etc. to reach Pico4-8-16 when capturing IR.
Trouble when IR-capture.	Variable strength or direction of transmitted IR. IR-signal not possible to copy.	Try different angles and distances (up to 10+ cm may be used) when doing IR-capture. Some IR-signals is not possible to copy.
When IR-capture: transmitter exits IR-capture when first IR-signal is sent.	Repetitive IR-signals from transmitter.	Avoid press <u>& hold</u> on key when IR-capture is done.

There are no fuses inside IR-transmitter.

7.1 If errors occur

If there are any faults that one not can find the answer for, please contact with your supplier for the equipment. Note details and phone number in table below.

Name	Telephone number
Company:	
Technician:	

8 Technical information

The Pico4, Pico8 and Pico16 IR-transmitters consists of a printed circuit board with a keyboard. These are placed inside a moulded plastic box. Serial number is noted on the printed circuit board, product name is engraved in the front cover. Battery is placed in the bottom part of the box. It is rechargeable and do not need to be replaced.

8.1 Battery

Please refer to chapter 2.7 and 2.7.1.

8.2 External switch

An external switch may be connected to the transmitter with a 3,5 mm mono jack. Switching is normally open, i.e. “N.O”. An external switch may be used in scanning mode. More information to be found in chapter 2.6.

8.3 Recycling

Electrical products shall be disposed according to national and international rules.

When a Picomed IR-transmitter is meant to be used by another person, the following should be done:

- ☐ Visual control of the transmitter: physical damages, cracks etc. in box or keypad?
- ☐ Cleaning.
- ☐ Functional check.
- ☐ Program an IR-receiver and the transmitter with the same IR-signals and test transmission.
- ☐ Connect an external switch and check scanning function.
- ☐ Execute IR-capture and check functionality.
- ☐ Recharge battery.

If transmitter shall be stored for a long time, it is advised to recharge battery to 50%: there is no indication for this so we advise to recharge it for approximately 1 hour before product is placed in stock. It is recommended to recharge it once a year.

8.4 Technical data


Type	Pico4, Pico8 and Pico16 IR-transmitter. Transmit built-in and copied IR-signals from other IR-transmitters. 2 pages.
Operation	Ordinary: 4/8/18 keys on keypad. Scanning: 3,5 mm mono jack for external switch, N.O. Light indicators in front. Programming: on external unit with screen via USB cable. Programming: some functions via keypad.
Number of channels	Pico4: 6. Pico8: 14. Pico16: 33.
Battery	Rechargeable, 3,7 VDC LiPo.
Battery life	Typically 3 months when in “ordinary use”.
Battery charging	4-5 hours. «USB» charger 5 V@0,5A, USB micro B.
Current consumption	<ul style="list-style-type: none"> • < 1,0 μA in standby. • $I_{RMS} \sim 25$ mA when IR-transmit.
Energy transmission	Infrared light, $\lambda = 940-950$ nm.
Operating temperature	-25 to +40° C.
Dimensions, weight	Pico4 and Pico8: W54-L120-H21 [mm], weight 80 gram.

Size of keys	Pico 16: W71-L150-H21 [mm], weight 140 gram.
Material in casing	Pico4 and Pico8 Ø = 14 mm, Pico16 Ø = 13 mm.
Transport	Plastic. Polycarbonate PC/TPU.
Development and manufacturing	To be packed in a way so no keys are pressed during transport.
	Norway, Picomed as, N-4993 SUNDEBRU.

8.5 CE-markings

IR transmitters is a subject to legislation in the EU and the EEA Agreement. They satisfy the requirements in “*Council Directive of 14 June 1993 concerning medical devices*”.

The product is in accordance with applicable EU/EEA directives and their related Norwegian regulations and they are CE marked. The corresponding Declaration of Conformity is in English version and is available upon request to manufacturer.



9 Appendix – programming map

